

Product Name :
Centrifuge, Digital RPM and Timer, Economy

Product Code :
ANALYTICAL SUK-AS03LAB050



Description :

Centrifuge, Digital RPM and Timer, Economy

Technical Specification :

Table best rotator with vigorous development and acrylic cover. With step less controller, Digital RPM meter and 0-99 minutes computerized clock and marker. Speed 0 to 4200 RPM. It can be fitted with various sort of rotor heads according to end-client prerequisites. Supply 230 volts/50Hz



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265