

**Product Name :**  
Centrifuge Machine Digital, 5200 r.p.m

**Product Code :**  
ANALYTICAL SUK-AS03LAB052



**Description :**

Centrifuge Machine Digital, 5200 r.p.m

**Technical Specification :**

Broadly useful Laboratory Centrifuge Machine, Microprocessor based (LCD Display). Overwhelming obligation C.R.C development completed with powder covering. Seat Top with Digital speedometer.digital clock of range 0-59 minutes and variable speed regulator.Maximum speed 5200 rpm(without stack). Fitted with light weight twofold walled top made of ABS plastic moulding.Fitted with very much adjusted all inclusive motor,Digital Speed Indicator, Programmable speed regulator,Lid bolt for safe drive



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265