

Product Name :
Micro Centrifuge, Micro Processor Based, 16000 rpm

Product Code :
ANALYTICAL SUK-AS03LAB053



Description :

Micro Centrifuge, Micro Processor Based, 16000 rpm

Technical Specification :

Rapid Micro Centrifuge Max rpm 16000, Heavy obligation plastic infusion formed body, Stainless Steel assurance bowl, Fitted with well adjust widespread engine having long existence with imported carbons, Microprocessor based Digital Timer of range 0-15 min. with LCD show, Preset speed controls through microchip from 2000 to 16000rpm, Digital show indicates Set Time/Run Time/rpm.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265