

**Product Name :**  
Micro Centrifuge, Micro Processor Based, 16000 rpm

**Product Code :**  
ANALYTICAL SUK-AS03LAB053



**Description :**

Micro Centrifuge, Micro Processor Based, 16000 rpm

**Technical Specification :**

Rapid Micro Centrifuge Max rpm 16000, Heavy obligation plastic infusion formed body, Stainless Steel assurance bowl, Fitted with well adjust widespread engine having long existence with imported carbons, Microprocessor based Digital Timer of range 0-15 min. with LCD show, Preset speed controls through microchip from 2000 to 16000rpm, Digital show indicates Set Time/Run Time/rpm.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265