

Product Name :
Refrigerated Micro Centrifuge, Digital, 16000 rpm

Product Code :
ANALYTICAL SUK-AS03LAB055



Description :

Refrigerated Micro Centrifuge, Digital, 16000 rpm

Technical Specification :

Table best refrigerated smaller scale axis machine. The temp. from room temp. upto short 8°C is controlled by PID electronic advanced temp. controller cum pointer. Microchip based Digital Timer of range 0-15 min. with LCD show. Preset speed controls through chip from 2000 to 16000 r.p.m . Advanced show demonstrates Set Time/Run Time/R.P.M. S.S. Bowl Size: 9" dia



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265