

**Product Name :**  
Table Top Centrifuge Machine High Speed - 20000 r.p.m.

**Product Code :**  
ANALYTICAL SUK-AS03LAB056



**Description :**

Table Top Centrifuge Machine High Speed - 20000 r.p.m.

**Technical Specification :**

Seat Top rapid Lab Centrifuge with advanced speedometer, computerized clock in the scope of 0-59 minutes, LED show and speed controller. Most extreme speed 20000 rpm when rotor head of 1.5ml of 12 tubes is fitted. Overwhelming obligation CRC development outfitted with powder covering. Fitted with all around adjusted widespread engine having long existence with imported carbons. Supply finished with clean cover, direction manual, rope and fitting to chip away at 220 volts 50 Hz A.C. This model is provided without rotor which can be chosen from beneath said alternatives according to client particular necessity.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

