

**Product Name :**  
Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m

**Product Code :**  
ANALYTICAL SUK-AS03LAB057



**Description :**

Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m

**Technical Specification :**

Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265