

Product Name :
Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m

Product Code :
ANALYTICAL SUK-AS03LAB057



Description :

Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m

Technical Specification :

Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265