

Product Name : Refrigerated Centrifuge, High Speed (Brushless Motors, Microprocessor based), 20000 rpm	Product Code : ANALYTICAL SUK-AS03LAB058
--	--



Description :

Refrigerated Centrifuge, High Speed (Brushless Motors, Microprocessor based), 20000 rpm

Technical Specification :

Special Features

Table top compact.

Brushless self balancing motor without carbons.

Ultra high speed.

Variety of rotors for different applications.

Elegant design.

S.S bowl.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265