

**Product Name :**  
Refrigerated Centrifuge Universal, 16000 rpm

**Product Code :**  
ANALYTICAL SUK-AS03LAB059



**Description :**

Refrigerated Centrifuge Universal, 16000 rpm

**Technical Specification :**

**Special Features**

Table Top centrifuge compact and elegant design.

Fitted with vacuum feet.

CFC free refrigerant.

Centrifuge bowl made of stainless steel.

Stable with no vibration.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265