

Product Name :
Refrigerated Centrifuge Universal, 16000 rpm

Product Code :
ANALYTICAL SUK-AS03LAB059



Description :

Refrigerated Centrifuge Universal, 16000 rpm

Technical Specification :

Special Features

Table Top centrifuge compact and elegant design.

Fitted with vacuum feet.

CFC free refrigerant.

Centrifuge bowl made of stainless steel.

Stable with no vibration.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265