

**Product Name :**  
Micro Centrifuge, High Speed, Brushless, 20000 rpm

**Product Code :**  
ANALYTICAL SUK-AS03LAB061



**Description :**

Micro Centrifuge, High Speed, Brushless, 20000 rpm

**Technical Specification :**

Rapid Micro Centrifuge, Max 20000 rpm, Heavy obligation C.R.C. development, Powder covered complete, Fitted with brushless programmed adjusted engine, Digital Timer of range 0-99 minutes, Digital Speed Indicator, Digital R.C.F. Marker, Programmable Speed controller, Electronic Safety top bolt, Imbalance identifier.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265