

**Product Name :**  
Hot Plate, Round Type, Thermostatic

**Product Code :**  
ANALYTICAL SUK-AS03LAB063



**Description :**

Hot Plate, Round Type, Thermostatic

**Technical Specification :**

A round cast press plate is settled on a thick mellow steel body. Body is covered with layer of paint. This give an exceedingly shiek and warmth safe surface for enduring appealing completion. Temperature is controlled by a vitality controller with constant kind warmth control. Rating 1.5 KW, 230 volts AC.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265