

**Product Name :**  
Hot Plate, Round Type, Digital

**Product Code :**  
ANALYTICAL SUK-AS03LAB064



**Description :**

Hot Plate, Round Type, Digital

**Technical Specification :**

A round cast press plate is settled on a thick gentle steel body. Body is secured with layer of paint. This give a profoundly shriek and warmth safe surface for dependable appealing completion. Temperature is controlled by a computerized temperature controller and pointer. Rating 1.5 KW, 230 volts AC



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265