

Product Name :
Hot Plate, Round Type, Digital

Product Code :
ANALYTICAL SUK-AS03LAB064



Description :

Hot Plate, Round Type, Digital

Technical Specification :

A round cast press plate is settled on a thick gentle steel body. Body is secured with layer of paint. This give a profoundly shriek and warmth safe surface for dependable appealing completion. Temperature is controlled by a computerized temperature controller and pointer. Rating 1.5 KW, 230 volts AC



Astrascientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265