

**Product Name :**  
Magnetic Stirrer without Hot Plate

**Product Code :**  
ANALYTICAL SUK-AS03LAB068



**Description :**

Magnetic Stirrer without Hot Plate

**Technical Specification :**

Maximum speed 1200 RPM. Suitable for stirring all volumes, supplied with a teflon coated magnetic stirring bar.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265