

**Product Name :**  
Magnetic Stirrer with Hot Plate, Digital RPM

**Product Code :**  
ANALYTICAL SUK-AS03LAB070



**Description :**

Magnetic Stirrer with Hot Plate, Digital RPM

**Technical Specification :**

Reasonable for blending all volumes. Most extreme speed 1200 rpm with variable speed control, provided with a teflon covered attractive mixing bar. Warmer of 300 watts with temperature controller. With Digital RPM marker.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265