

**Product Name :**  
Laboratory Stirrer, Economy

**Product Code :**  
ANALYTICAL SUK-AS03LAB073



**Description :**

Laboratory Stirrer, Economy

**Technical Specification :**

With clamor less widespread engine, greatest speed 4000 rpm with electronic variable speed controller between 200 to 4000 rpm. The stirrer shaft is of stainless steel with 3 bladed propeller having 37 mm. breadth, shaft and stirrer coupled by methods for straightforward shrubbery coupling with grub screw. Stirrer is mounted on a customizable overwhelming stand.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265