

Product Name :
Laboratory Stirrer, High Torque, Geared Motor

Product Code :
ANALYTICAL SUK-AS03LAB074



Description :

Laboratory Stirrer, High Torque, Geared Motor

Technical Specification :

Safe and a versatile unit. A very high torque is developed at a relatively low operating speed. It is suitable for affectively stirring high viscosity fluids and pastes.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265