

Product Name :
Water Bath, Double Walled, Digital Controller

Product Code :
ANALYTICAL SUK-AS03LAB079



Description :

Water Bath, Double Walled, Digital Controller

Technical Specification :

Water Bath, twofold walled, advanced controller is reasonable for vanishing fluids like liquor or for softening solids, for example, agar. Concentric rings can be expelled to accomodate compartments of different sizes. Development: Double walled chambers made of stainless steel and external divider made of thick gentle steel completed with epoxy based powder covered paint. The twofold walled chamber is loaded with high-review protection to avoid warm misfortune. This is given 1500 watts gadget controlled by electronic computerized temperature controller cum pointer from surrounding to 110°C within+/- 1°C, two marker lights and one On/off Rotary Switch

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com