

Product Name :
Water Bath, Double Walled, Stainless Steel, Digital Controller

Product Code :
ANALYTICAL SUK-AS03LAB080



Description :

Water Bath, Double Walled, Stainless Steel, Digital Controller

Technical Specification :

Water Bath, twofold walled, stainless steel computerized controller is appropriate for dissipating fluids like liquor or for softening solids, for example, agar. Concentric rings can be expelled to accommodate compartments of different sizes. Development: Double walled chambers made of stainless steel and external divider made of 304 review Stainless Steel completed with epoxy based powder covered paint. The twofold walled chamber is loaded with high-review protection to avoid warm misfortune. This is furnished with 1500 watts gadget controlled by electronic computerized temperature controller cum marker from surrounding to 110°C within+/- 1°C, two pointer lights and one On/off Rotary Switch.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265