

Product Name :
Serological Water Bath, Digital Temperature Controller, All Stainless Steel

Product Code :
ANALYTICAL SUK-AS03LAB083



Description :

Serological Water Bath, Digital Temperature Controller, All Stainless Steel

Technical Specification :

Serological water bath is suitable for clinical procedure, incubation, inactivation & pharmaceutical biomedical procedure.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265