

Product Name :
Cylindrical Oil Bath (High Temperature)

Product Code :
ANALYTICAL SUK-AS03LAB089



Description :

Cylindrical Oil Bath (High Temperature)

Technical Specification :

Twofold Walled. With temperature extend 50 to 300 degree celsius
Controlled by Thermostat/Digital Controller. Inward chamber made of
stainless steel. Accessible both with indoor regulator and advanced
temperature controller cum marker.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265