

**Product Name :**  
Ultra Sonic Bath (Sonicator, Ultra Sound)

**Product Code :**  
ANALYTICAL SUK-AS03LAB092



**Description :**

Ultra Sonic Bath (Sonicator, Ultra Sound)

**Technical Specification :**

It is the sound which person can't here Ultrasound utilizes frequencies from and above 18KHz (18,000 cycles for each stable). High recurrence electrical vitality is changed over into ultrasound waves by methods for ultrasonic transducers, which are fortified on the base of S.S. cleaning tank. These high frequency sound waves make in the fluid incalculable minuscule vacuum bubbles which quickly extend and crumple.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265