

**Product Name :**  
Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel

**Product Code :**  
ANALYTICAL SUK-AS03LAB095



**Description :**

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel

**Technical Specification :**

Available Sizes:

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 2'x2'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 4'x 2'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 4'x 3'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 4'x4'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 6'x4'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 8'x8'

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

