

Product Name :
Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel

Product Code :
ANALYTICAL SUK-AS03LAB095



Description :

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel

Technical Specification :

Available Sizes:

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 2'x2'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 4'x 2'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 4'x 3'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 4'x4'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 6'x4'

Laminar Air Flow, Vertical, Suspended Ceiling Mounted, Stainless Steel, Working Area 8'x8'

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

