

**Product Name :**  
Rotary / Orbital Shaker

**Product Code :**  
ANALYTICAL SUK-AS03LAB102



**Description :**

Rotary / Orbital Shaker

**Technical Specification :**

Offbeat drive bestows uniform round movement to the stage in a level plane. Sturdy all steel development, completed in blue stoved finish. Finish with engine to take a shot at 220 volts AC. It has control for variable speed. Accessible in various designs with advanced clock and computerized RPM meter on request.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265