

**Product Name :**  
VDRL Shaker

**Product Code :**  
ANALYTICAL SUK-AS03LAB103



**Description :**

VDRL Shaker

**Technical Specification :**

It is extremely helpful for V.D.R.L. gathering and writing and in addition other agglutination tests. It is likewise utilized for blending of arrangement in little carafe, measuring glasses, bottles or comparative compartments. It has a clock consolidated in the electrical circuit which is being cut off naturally when the setting is finished. Stage 30x30 cm. Likewise accessible with advanced clock.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265