

**Product Name :**  
Flask Shaker Wrist Action Type

**Product Code :**  
ANALYTICAL SUK-AS03LAB104



**Description :**

Flask Shaker Wrist Action Type

**Technical Specification :**

This model is made from heavy cast iron & finished with a thick layer of epoxy based polyester paint.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265