

**Product Name :**  
Willy Grinder Or Lab Willy Mill (Arthur Thomas Type)

**Product Code :**  
ANALYTICAL SUK-AS03LAB105



**Description :**

Willy Grinder Or Lab Willy Mill (Arthur Thomas Type)

**Technical Specification :**

Accessible in different chamber sizes with cover. Three conveyance strainers 20, 40 and 60 work individually, each with plastic collector. Takes a shot at 220 volts AC. A removeable glass plate shapes the substance of the chamber and allows perception amid activity.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265