

**Product Name :**  
Spectrophotometer, Digital, Double Display

**Product Code :**  
ANALYTICAL SUK-AS03LAB108



**Description :**

Spectrophotometer, Digital, Double Display

**Technical Specification :**

Range: 340 to 960 nm

Resolution: 0.001 Abs

Output: 0 to 100%, Abs; 0 to 1.999, Conc; 0 to 1999

Detector: Photodiode

Light Source: Vertical tungsten halogen lamp

Monochromator: Integrating

Path Length:

Spectral Bandwidth:

1 nm

Resolution:

Wave Length

Wave Length

Accuracy:

Repeatability:

Display LED display

Operating Temp.:

Power 10% AC 50 Hz

Dimensions: 165 mm

Weight (approx.)



Astra Scientific

## Astra Scientific

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265