


Product Name :
Spectrophotometer, Microprocessor Based, VIS, Single Beam

Product Code :
ANALYTICAL SUK-AS03LAB110



Description :

Spectrophotometer, Microprocessor Based, VIS, Single Beam

Technical Specification :

- Range 320 to 1020 nm.
- Absorbance and K factor measurement.
- 3 ½ bright LED display.
- Unlimited sample data storage facility.
- Printer attachment interface via RS 232C.
- Accommodate 10 mm & 50 mm path length cuvettes.
- Supplied with 4 glass cuvettes of 10 mm path length and software



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265