

Product Name : Spectrophotometer, Microprocessor Based, UV-VIS Single Beam, Automatic	Product Code : ANALYTICAL SUK-AS03LAB112
---	--



Description :

Spectrophotometer, Microprocessor Based, UV-VIS Single Beam, Automatic

Technical Specification :

~~Single Beam~~ Grating 1200 lines/mm

~~Wavelength~~ Range:

~~Band~~ Width:

~~Wavelength~~ Accuracy:

~~Wavelength~~ Repeatability:

~~Wavelength~~ Setting:

~~Photometric~~ Accuracy:

+0.3 %T

Repeatability:

Photometric

Thermistor Accuracy 0-9999C

Range:

Stability @ 500nm

Stray Light @ 220nm, 360nm

USB Output Port:

Parallel Port

Display: Dots LCD

Deuterium Lamp & Tungsten Halogen Lamp

Silicon Photodiode

Power Requirement:



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265