

Product Name :
Spectrophotometer, Atomic Absorbtion

Product Code :
ANALYTICAL SUK-AS03LAB116



Description :

Spectrophotometer, Atomic Absorbtion

Technical Specification :

Wavelength range:

Wavelength accuracy:
± 0.25 nm

Resolution:
spectral lines of Mn at 279.5 nm and 279.8 nm can be separated with the spectral bandwidth of 0.2 nm an valley/peak energy ratio less than 30%.

Baseline stability:

Back ground correction:

- The D-2 lamp background correction capability at 1A is better than 30 times.
- The S-H background correction capability at 1.8 A is better than 30

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

