

Product Name : Spectrophotometer, Atomic Absorbtion	Product Code : ANALYTICAL SUK-AS03LAB116
---	--



Description :
Spectrophotometer, Atomic Absorbtion
Technical Specification :
Wavelength range: 190-900 nm
Wavelength accuracy: ± 0.25 nm
Resolution: spectral lines of Mn at 279.5 nm and 279.8 nm can be separated with the spectral bandwidth of 0.2 nm an valley/peak energy ratio less than 30%.
Baseline stability: 0.04% / 30 min
Back ground correction: <ul style="list-style-type: none">• The D-2 lamp background correction capability at 1A is better than 30 times.• The S-H background correction capability at 1.8 A is better than 30

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

