

Product Name :
PH Meter, Digital, A.T.C. Waterproof (Pocket Model)

Product Code :
ANALYTICAL SUK-AS03LAB118



Description :

PH Meter, Digital, A.T.C. Waterproof (Pocket Model)

Technical Specification :

Range pH: 0.0 – 14.0 pH

Range Temp.: 0°C – 55°C

Resolution pH: 0.01 pH

Resolution Temp.: 0.1°C / 0.1°F

Accuracy pH: ± 0.1 pH

Accuracy Temp.: ± 1°C

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific