

Product Name :
PH Meter, Digital, Table Model, Automatic Temperature Compensation

Product Code :
ANALYTICAL SUK-AS03LAB122



Description :

PH Meter, Digital, Table Model, Automatic Temperature Compensation

Technical Specification :

Display: 3½ digit LCD

Range: pH : 0 to 14.00

mV : 0 to \pm 1999

Resolution: pH : 0.01

mV : 1

Accuracy: pH : \pm 0.01

mV : 1

Temperature: Auto 0 to 100°

Compensation: Automatic 0 to 100°

Slope Correction: 80 to 120%

Power: 230V \pm 10%

Dimensions: 76 x 275 x 175 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265