

**Product Name :**  
PH Meter, Digital, Table Model, Automatic Temperature Compensation

**Product Code :**  
ANALYTICAL SUK-AS03LAB122



**Description :**

PH Meter, Digital, Table Model, Automatic Temperature Compensation

**Technical Specification :**

Display: 3½ digit LCD

Range: pH : 0 to 14.00

mV : 0 to ± 1999

Resolution: pH : 0.01

mV : 1

Accuracy: pH : ±0.01

mV : 1

Temperature: Auto 0 to 100°

Compensation: Automatic 0 to 100°

Slope Correction: 80 to 120%

Power: 230V  $\pm$ 10%

Dimensions: 76 x 275 x 175 mm



## **Astra Scientific**

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265