

Product Name :
PH, Conductivity and Temperature Meter, Digital, Table Model

Product Code :
ANALYTICAL SUK-AS03LAB123



Description :

PH, Conductivity and Temperature Meter, Digital, Table Model

Technical Specification :

General

Power: 220 V AC \pm 10%, 50Hz

Dimensions: 75 x 175 x 275 mm (App.)

Weight: 2.5 Kg. (App.)

Temp. Compensation: 0 to 50°C adjustable from temperature.

Parameters

pH Range: 0 - 14 pH

Resolution: 0.01 pH

Accuracy: \pm 0.01 pH \pm 1 Digit

Temp. Compensation Manual: 0 to 100°C

Temperature

Range: 0 - 100°C

Resolution: 0.1°C

Accuracy: $\pm 0.2\% \pm 1$ Digit

Sensor RTD: (PT - 100)

Conductivity

Range: 0 to 20 mS/cm in two Ranges

Accuracy: $\pm 0.5\% \pm 1$ Digit

Resolution: 1 μ S/cm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265