

**Product Name :**  
PH, Conductivity and Temperature Meter, Digital, Table Model

**Product Code :**  
ANALYTICAL SUK-AS03LAB123



**Description :**

PH, Conductivity and Temperature Meter, Digital, Table Model

**Technical Specification :**

**General**

Power: 220 V AC  $\pm$  10%, 50Hz

Dimensions: 75 x 175 x 275 mm (App.)

Weight: 2.5 Kg. (App.)

Temp. Compensation: 0 to 50°C adjustable from temperature.

**Parameters**

pH Range: 0 - 14 pH

Resolution: 0.01 pH

Accuracy:  $\pm$  0.01 pH  $\pm$  1 Digit

Temp. Compensation Manual: 0 to 100°C

#### Temperature

Range: 0 - 100°C

Resolution: 0.1°C

Accuracy:  $\pm 0.2\% \pm 1$  Digit

Sensor RTD: (PT - 100)

#### Conductivity

Range: 0 to 20 mS/cm in two Ranges

Accuracy:  $\pm 0.5\% \pm 1$  Digit

Resolution: 1 $\mu$ S/cm



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265