

Product Name :
PH/Temperature and mV Meter, Microcontroller Based

Product Code :
ANALYTICAL SUK-AS03LAB125



Description :

PH/Temperature and mV Meter, Microcontroller Based

Technical Specification :

Range: 4 pH

Resolution:

Accuracy: ± 1 digit

80% to 120% (default 100%).

Temp: 0 to 100°C (with temp. probe/Manual with Keys).

Auto Buffer: 0 to 9.20 pH (AUTOMATIC)

Auto & MANUAL, 3 Point/2 Point (Programmable) 2- POINT, 3- POINT

pH error, ORP error, CAL error & T- error

Temperature

Range: 0 +200 °C

Resolution:

Accuracy: of range + 1 digit

Probe Type

Display: touch keys with audible sound

Storage

Real time display with Date

Ports: RS 232 & Parallel

Printers: Date Wise (Programmable)

Digital Output: & Printer Centronics (Programmable)

Error Indication:



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265

