

Product Name : PH/Temperature and mV Meter, Microcontroller Based 2	Product Code : ANALYTICAL SUK-AS03LAB128
-------------------------------------------------------------------------------	----------------------------------------------------



Description :

PH/Temperature and mV Meter, Microcontroller Based 2

Technical Specification :

Technical Specification

- Range: 0-14 pH
- Resolution: 0.01 pH
- Accuracy: ± 0.01 pH
- 80% to 120% of range
- Temperature Compensation: temp. probe/Manual with Keys.
- Auto Buffer Replenish: 0-20 pH (AUTOMATIC)
- Calibration: MANUAL, 3 Point/2 Point (Programmable) 2- POINT, 3- POINT
- Temperature Range: 0 +200 °C
- Resolution: 0.1 °C
- Accuracy: ± 0.1 °C
- Probe Type: Pt100 probe
- Data Entry: touch keys with audible sound
- 5000 Samples
- Real time display with Date
- RS 232 & Parallel
- Print Data & Date Wise (Programmable)
- Digital RS 232 & Printer Centronics (Programmable)
- Error Indication: Error, CAL error & T- error



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265