

Product Name :
Computerized Bench Top Cold Storage Trainer

Product Code :
HEAT SUK-AS03LAB006



Description :

Computerized Bench Top Cold Storage Trainer

Technical Specification :

Computerized Bench Top Cold Storage Trainer: refrigeration lab equipment: The trainer features a large refrigeration chamber. Two evaporators allow for an investigation of the effect of different evaporator sizes on the cold storage room climate and the icing. An electric defrost heater and hot gas defrosting



Astrascientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265