

Product Name :
Computerized Cascade Refrigeration Trainer

Product Code :
HEAT SUK-AS03LAB007



Description :

Computerized Cascade Refrigeration Trainer

Technical Specification :

Computerized Cascade Refrigeration Trainer: air conditioning trainer: The cascade Refrigeration systems with two-stage compression are used for the generation of particularly low temperatures. At very low temperatures large pressure differences are required between the evaporator and condenser. In a compressor



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265