

Product Name :
Computerized Cooling Tower Test Bench

Product Code :
HEAT SUK-AS03LAB008



Description :

Computerized Cooling Tower Test Bench

Technical Specification :

Computerized Cooling Tower Test Bench: air conditioning trainer: The examines the main components and principle of a wet cooling tower with forced ventilation. Water is heated in a tank and transported by a pump to an atomizer. The atomizer sprays the water to be cooled over the wet deck surface.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265