

**Product Name :**  
Vapour Compression Refrigeration Trainer

**Product Code :**  
HEAT SUK-AS03LAB013



**Description :**

Vapour Compression Refrigeration Trainer

**Technical Specification :**

Vapour Compression Refrigeration Trainer: The experimental unit represents a typical refrigeration circuit consisting of a hermetic piston compressor, condenser, expansion valve and evaporator. The evaporator and condenser are transparent to provide good monitoring of the phase transition process during evaporation



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265