

Product Name :
Vapour Jet Refrigeration Trainer

Product Code :
HEAT SUK-AS03LAB014



Description :

Vapour Jet Refrigeration Trainer

Technical Specification :

Vapour Jet Refrigeration Trainer: The Unlike standard compression refrigeration systems, vapour jet refrigerating machines do not have a mechanical but a vapour jet compressor. This makes it possible to use different heat sources for cold production. Such sources could e.g. be solar energy or process waste heat



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265