

<b>Product Name :</b> Cut Section Model Of Open Hermetic Compressor	<b>Product Code :</b> HEAT SUK-AS03LAB015
	
<b>Description :</b>	
Cut Section Model Of Open Hermetic Compressor	
<b>Technical Specification :</b>	
Cut Section Model Of Open Hermetic Compressor: Refrigeration and Air conditioning lab equipment: Cut Section Model Of Open Hermetic Compressor model will be made out of original parts. The model will be suitably sectioned to demonstrate the internal construction details showing the minute information	



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265