

Product Name :
Taper Gauge

Product Code :
MEASUREMENT SUK-
AS03LAB024



Description :

Taper Gauge

Technical Specification :

Taper Gauge: analytical instrument: The taper plug gage is an internal gage in the shape of the frustum of a cone. It is used to measure internal tapers. It has wide application in design engineering, mechanical



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265