

Product Name :
Computerized Gear Pump Test Rig

Product Code :
TOURBO SUK-AS03LAB002



Description :

Computerized Gear Pump Test Rig

Technical Specification :

Computerized Gear Pump Test Rig: The designed and developed gear pump is the most widely used of the positive action rotary pumps. Two gear wheels operate inside a casing. One is driven while the other rotates in mesh with it. The liquid is carried around in the space between consecutive teeth and then ejected as the teeth mesh. The pump has no valves. It is a positive displacement



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265